



KEBBY INDUSTRIES

4075 KILBURN AVENUE
ROCKFORD, IL 61101 U.S.A.

Following is a list of the materials in Kebby Hand Crimpers and Decappers:

INTERNAL PARTS:

Stainless Steel
Carbon Steel
Delrin
Lubricant: Food grade grease

MATERIALS CONTACTING SEAL DURING CRIMP PROCESS:

Stainless Steel
Nickel-Plated Carbonized Steel

RECOMMENDED CLEANING:

Wipe with a clean cloth or mild cleanser

STORAGE:

Store in box provided in dry area