



KEBBY INDUSTRIES

4075 KILBURN AVENUE
ROCKFORD, IL 61101 U.S.A.

Materials for parts in Kebby products:

Power Crimper:

- Exterior parts are Electroless nickel-plated, aluminum, and stainless steel.
- Interior parts are aluminum, stainless steel, Delrin, Buna-N, Teflon, and urethane.

Power Crimper and Decrimper Head:

- Aluminum
- Electroless Nickel Plated parts
- Stainless Steel
- Steel
- Delrin

Hand Crimper and Decrimper, Painted:

- **INTERNAL PARTS:**
Stainless Steel
Carbon Steel
Delrin
Lubricant: Food grade grease
- **MATERIALS THAT HAVE CONTACT WITH SEALS:**
Stainless Steel
Nickel-Plated Carbonized Steel

s:klo/crimpers decappers/materials for Kebby items